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Flaminal® Hydro

Properties

Flaminal® Hydro is a new class of wound healing agent, called enzyme alginate. It consists of an antibacterial enzyme system embedded in hydrated alginates. Flaminal® Hydro supports the healing of a wound in a number of ways:

1. Flaminal® Hydro keeps the wound clean.
2. Flaminal® Hydro keeps the wound moist.
3. Flaminal® Hydro offers antimicrobial protection.
4. Flaminal® Hydro is safe for the skin and wound tissue.
5. Flaminal® Hydro reduces wound odour caused by bacteria.
6. Flaminal® Hydro helps to reduce excessive protease activity (protease = enzyme that breaks down proteins).
7. Flaminal® Hydro reduces bacteria released from a biofilm.
8. Flaminal® Hydro protects the wound border.


These are the prime properties for ensuring rapid wound healing. Flaminal® Hydro is indicated for slightly to moderately exuding wounds such as:

- leg ulcers
- diabetic ulcers
- 2nd degree burns (deep, superficial)
- wounds from radiotherapy
- oncology wounds
- complex grazes
- wounds from dermatosurgery

How to use

Applying or renewing Flaminal® Hydro is easy and painless. Proceed as follows:

1. Rinsing and cleansing the wound and surrounding area is recommended prior to application of Flaminal® Hydro.
2. Cover the entire wound bed with a sufficiently thick layer (e.g. 4-5 mm) of Flaminal® Hydro
3. Cover Flaminal® Hydro with a dressing. The choice of dressing should be based on the type of lesion:
 - Dry wound: transparent film (polyurethane) or paraffin gauze fixed by a non-adhesive bandage
 - Slightly exuding: non-adherent dressing fixed by a non-adhesive bandage or by a hypoallergenic, adhesive tape
 - Moderately exuding: absorbent, non-adherent compress fixed by a non-adhesive bandage or a hypoallergenic, adhesive tape

Good to know

- Check the dressing daily. If it is leaking or if there is insufficient gel, Flaminal® Hydro should be changed and a new dressing applied.
- During Flaminal® Hydro treatment, dry whitish alginate flakes may appear on the wound border. They should not be removed. The flakes will prevent wound border maceration. Maceration slows the healing process of the wound.
- If dry whitish alginate flakes appear in the wound, cover Flaminal® Hydro as you would with a dry wound. The flakes in the wound will disappear when wound moisture balance is restored.
- Should Flaminal® Hydro liquefy too quickly, the wound is too exuding for Flaminal® Hydro. Use Flaminal® Forte.
- The dressing with Flaminal® Hydro can remain in place as long as the gel structure is intact (1 to 4 days, depending on the amount of exudate).
- In the initial phase of the treatment the wound may seem to become larger. This is a normal part of the healing process in which dead tissue is eliminated.
- The colour of the alginates may change slightly due to their natural origin. Colour changes have no effect on the quality of Flaminal® Hydro.

Special wounds

- Leg ulcers: adhesives often cause allergic reactions, particularly in patients who suffer from leg ulcers. Use a non-adhesive bandage.
- Undermined wounds: put sterile gauze tampons in the cavities as a support and apply a sufficient amount of Flaminal® Hydro in the cavities between the tampons.
- Fistulas and small deep cavities: inject Flaminal® Hydro carefully into the cavity using a dedicated cannula.
- Exuding wounds: apply a thicker layer of Flaminal® Hydro or use Flaminal® Forte.
- Heel ulcers: apply Flaminal® Hydro and cover with a V-shaped gauze compress.

Precautions

Some wounds may indicate a complementary treatment. Consult a physician or wound specialist. Flaminal® Hydro can be used on infected wounds under medical supervision. Alginates have a faint odour. This odour has no influence on the healing process. Should there be a sudden distinct odour, consult a physician or wound specialist. Flaminal® Hydro is not indicated in patients with a known allergy to one of the components (see under "Composition").

Flaminal® Hydro cannot be used on the eyelids or in the eye. Should it come into contact with an eye, rinse the eye thoroughly with running water and consult a physician.

Storage

Keep at room temperature (below 25°C) in a dry place in the original pack. Do not use Flaminal® Hydro after the expiry date shown on the pack. The expiry date (USE BY: Year/Month) is the last day of the indicated month and year. Recap the tube or jar carefully immediately after use. Once opened, and if recapped carefully, a tube of Flaminal® Hydro can be stored and used until the expiry date on the tube and box. A jar of Flaminal® Hydro 500g however cannot be used more than one week after its first opening. We recommend that a pack of Flaminal® Hydro be used for one patient only (single patient).

Composition

Alginate (total alginate content 3.5% (g/g)); macrogol; hydroxypropyl-cellulose; enzyme system (glucose oxidase, lactoperoxidase, glucose, guaiacol); potassium sorbate; potassium iodide; buffer; purified water Ph.Eur.

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- UK: Flen Health Ltd. www.flenhealth.co.uk